



fitted with	1 x LED module - Energy efficiency category A+
Illuminancefield @ 100 cm	Ø 16 cm
illuminance @ 100 cm	60 000 lx
color rendering index (CRI)	> 95
colour temperatur	4500 K
Depth of illumination L1+L2	1.2 m
work equipment	electronic converter
mains lead	direct power supply
luminaire body	polycarbonate & acrylonitrile butadiene styrene copolymer (PC & ABS)
weight (net)	13.0 kg
Power consumption	30 W
usage	rotary switch on the luminaire head single-pole secondary, continuously dimmable spring (weight compensation)
balance of articulation joints	I
class of protection	I
Accessory	sterilizable handle
norms	EN 60601-1, EN 60601-2-41
special features	sterilisable grip, continuously dimmable, maintenance free
Fastening	wall mount
System of protection	IP 43 (Leuchtenkopf)
order number	D15175000 - pure white

light distribution

