



fitted with	1 x LED module - Energy efficiency category A+
Illuminancefield @ 50 cm	Ø 18 cm
illuminance @ 50 cm	40 000 lx
color rendering index (CRI)	> 93
colour temperatur	4400 K / 3800 K / 3300 K
work equipment	electronic converter
mains lead	1,8 m Euro
luminaire body	polycarbonate (PC)
weight (net)	1.4 kg
Power consumption	18 W
usage	keypad at luminaire head
balance of articulation joints	friction
class of protection	II
Accessory	roll support, rail fastening with removal protection, table clamp with removal protection, universal fastening with removal protection, extension arm with removal protection, wall bracket with removal protection
norms	EN 60601-1, EN 60601-2-41
special features	3 color shades, 5 dimming levels, maintenance free
Fastening	rail, table, wall
order number	D15473000 - pure white

light distribution

measuring distance 0.5 m

