



| | |
|---------------------------------------|---|
| fitted with | 1 x LED module - Energy efficiency category A+ |
| Illuminancefield @ 50 cm | Ø 18 cm |
| illuminance @ 50 cm | 40 000 lx |
| color rendering index (CRI) | > 93 |
| colour temperatur | 4400 K |
| work equipment | electronic converter |
| mains lead | 1,8 m Euro |
| luminaire body | polycarbonate (PC) |
| weight (net) | 1.4 kg |
| Power consumption | 13 W |
| usage | rocker switch (I/O) at the arm |
| balance of articulation joints | spring |
| class of protection | II |
| Accessory | roll support, rail fastening with removal protection, table clamp with removal protection, universal fastening with removal protection, extension arm with removal protection, wall bracket with removal protection |
| norms | EN 60601-1, EN 60601-2-41 |
| special features | switchable, maintenance free |
| Fastening | rail, table, wall |
| order number | D15472000 - pure white |



light distribution

measuring distance 0.5 m

