



fitted with	1 x LED module - Energy efficiency category A+
Illuminancefield @ 100 cm	Ø 17 cm
illuminance @ 100 cm	2 x 50000 lx
color rendering index (CRI)	> 95
colour temperatur	4300 K
Depth of illumination L1+L2	1,47 m
work equipment	electronic converter
connected load	230-240 V; 50-60 Hz
mains lead	direct power supply
luminaire body	polycarbonate & acrylonitrile butadiene styrene copolymer (PC & ABS)
weight (net)	22.0 kg
Power consumption	32 W
balance of articulation joints	spring
class of protection	I
Accessory	sterilizable handle, adapter for suspended ceiling
norms	EN 60601-1, EN 60601-2-41
special features	with grip, maintenance free
Fastening	ceiling mount
System of protection	IP 20
order number	D15612000 - pure white

light distribution

