



fitted with	1 x LED module - Energy efficiency category A+
Illuminancefield @ 100 cm	Ø 17 cm
illuminance @ 100 cm	50 000 lx
color rendering index (CRI)	> 95
colour temperatur	4300 K
Depth of illumination L1+L2	1,47 m
work equipment	electronic converter
mains lead	direct power supply
luminaire body	polycarbonate & acrylonitrile butadiene styrene copolymer (PC & ABS)
weight (net)	14.0 kg
Power consumption	16 W
usage	rocker switch (I/O) at the arm, 2-stages rocker switch (60/100 %)
balance of articulation joints	spring (weight compensation)
class of protection	I
Accessory	sterilizable handle
norms	EN 60601-1, EN 60601-2-41
special features	with grip, maintenance free
Fastening	wall mount
System of protection	IP 20
order number	D15613000 - pure white

light distribution

