



|                                       |  |
|---------------------------------------|--|
| <b>fitted with</b>                    | 1 x LED module - Energy efficiency category A+                       |
| <b>Illuminancefield @ 100 cm</b>      | Ø 17 cm  |
| <b>illuminance @ 100 cm</b>           | 50 000 lx  |
| <b>color rendering index (CRI)</b>    | > 95   |
| <b>colour temperatur</b>              | 4300 K   |
| <b>Depth of illumination L1+L2</b>    | 1,47 m   |
| <b>work equipment</b>                 | electronic converter   |
| <b>mains lead</b>                     | direct power supply  |
| <b>luminaire body</b>                 | polycarbonate & acrylonitrile butadiene styrene copolymer (PC & ABS) |
| <b>weight (net)</b>                   | 13.0 kg  |
| <b>Power consumption</b>              | 16 W   |
| <b>usage</b>                          | rocker switch (I/O) at the arm, 2-stages rocker switch (60/100 %)    |
| <b>balance of articulation joints</b> | spring (weight compensation)   |
| <b>class of protection</b>            | I  |
| <b>Accessory</b>                      | sterilizable handle, adapter for suspended ceiling                   |
| <b>norms</b>                          | EN 60601-1, EN 60601-2-41  |
| <b>special features</b>               | with grip, maintenance free  |
| <b>Fastening</b>                      | ceiling mount  |
| <b>System of protection</b>           | IP 20  |
| <b>order number</b>                   | D15611000 - pure white   |

### light distribution

