



fitted with	LED - Energy efficiency category A+
color rendering index (CRI)	> 95
colour temperatur	6500 K
work equipment	electronic converter
mains lead	1.8 m Euro
lens; dimensions	Ø 160 mm
lens; diopters	3,5
luminaire body	polyamide (PA)
weight (net)	2.1 kg
Power consumption	12 W
usage	keypad at luminaire head
balance of articulation joints	spring
class of protection	II
Accessory	roll support, rail fastening with removal protection, table clamp with removal protection, universal fastening with removal protection, extension arm with removal protection, wall bracket with removal protection
norms	EN 60601-1, EN 60601-2-41
special features	Rotary feedthrough (360° Turnover axis), scratch-proof plastic lens, continuously dimmable
Fastening	roller stand, rail, table, wall mount, pin (16 mm W)
order number	D15815000 - pure white



light distribution

