



fitted with	LED - Energy efficiency category A+> 95
color rendering index (CRI)	> 95
colour temperatur	6500 K
work equipment	electronic converter
mains lead	direct power supply
lens; dimensions	Ø 160 mm
lens; diopters	3,5
luminaire body	polyamide (PA)
weight (net)	9.1 kg
Power consumption	12 W
usage	keypad at luminaire head
balance of articulation joints	spring
class of protection	I
Accessory	adapter for suspended ceiling
norms	EN 60601-1, EN 60601-2-41
special features	Rotary feedthrough (360° Turnover axis), scratch-proof plastic lens, continuously dimmable
Fastening	ceiling mount, ceiling tube can be cut to length (40-116cm)
order number	D15817000 - pure white



light distribution

