



| | |
|---------------------------------------|---|
| fitted with | LED - Energy efficiency category A+ > 95 |
| color rendering index (CRI) | > 95 |
| colour temperatur | 6500 K |
| work equipment | electronic converter |
| mains lead | direct power supply |
| lens; dimensions | Ø 160 mm |
| lens; diopters | 3,5 |
| luminaire body | polyamide (PA) |
| weight (net) | 9.1 kg |
| Power consumption | 12 W |
| usage | keypad at luminaire head |
| balance of articulation joints | spring |
| class of protection | I |
| Accessory | adapter for suspended ceiling |
| norms | EN 60601-1, EN 60601-2-41 |
| special features | Rotary feedthrough (360° Turnover axis), scratch-proof plastic lens, continuously dimmable, switchable Wood-Light |
| Fastening | ceiling mount, ceiling tube can be cut to length (40-116cm) |
| order number | D15816000 - pure white |

light distribution

