



fitted with	LED - Energy efficiency category A+
Illuminancefield @ 50 cm	Ø 18.0 cm
illuminance @ 50 cm	50 000 lx
color rendering index (CRI)	> 93
colour temperatur	4400 K
work equipment	electronic converter
mains lead	1.8 m Euro
luminaire body	PMMA
weight (net)	1.1 kg
Power consumption	13 W
usage	rocker switch (I/O) at the arm
class of protection	II
Accessory	roll support, rail fastening with removal protection, table clamp with removal protection, extension arm with removal protection, wall bracket with removal protection
norms	EN 60601-1, EN 60601-2-41
special features	with grip
Fastening	roller stand, rail, table, wall, pin (16 mm W)
System of protection	IP 20
order number	D15732000 - pure white



light distribution

