



|                                    |  |
|------------------------------------|--|
| <b>fitted with</b>                 | LED - Energy efficiency category A+  |
| <b>Illuminancefield @ 50 cm</b>    | Ø 18.0 cm  |
| <b>illuminance @ 50 cm</b>         | 50 000 lx  |
| <b>color rendering index (CRI)</b> | > 93   |
| <b>colour temperatur</b>           | 4400 K   |
| <b>work equipment</b>              | electronic converter   |
| <b>mains lead</b>                  | 1.8 m US   |
| <b>luminaire body</b>              | PMMA   |
| <b>weight (net)</b>                | 1.1 kg   |
| <b>Power consumption</b>           | 13 W   |
| <b>usage</b>                       | rocker switch (I/O) at the arm   |
| <b>class of protection</b>         | II   |
| <b>Accessory</b>                   | roll support, rail fastening with removal protection, table clamp with removal protection, extension arm with removal protection, wall bracket with removal protection |
| <b>norms</b>                       | EN 60601-1, EN 60601-2-41  |
| <b>special features</b>            | with grip  |
| <b>Fastening</b>                   | roller stand, rail, table, wall, pin (16 mm W)   |
| <b>System of protection</b>        | IP 20  |
| <b>order number</b>                | D15757000 - pure white   |



## light distribution

